

## (7) ALPS処理水放出口北2km西0.5km

採取日	核種濃度 (Bq/L ※PuはmBq/L)									備考
	全β放射能	<sup>134</sup> Cs	<sup>137</sup> Cs	<sup>3</sup> H			<sup>90</sup> Sr	<sup>238</sup> Pu	<sup>239+240</sup> Pu	
				減圧蒸留法	迅速分析	電解濃縮法				
R4. 5. 19	0.02	ND	0.012	ND		0.10	0.0010	ND	ND	
R4. 8. 2	0.02	ND	0.005	ND		0.10	0.0007	ND	ND	
R4. 11. 8	0.03	ND	0.015	ND		0.06	0.0011	ND	ND	
R5. 2. 7	0.01	ND	0.007	ND		0.05	0.0010	ND	ND	
R5. 5. 10	0.01	ND	0.009	ND		ND	0.0009	ND	0.010	
R5. 8. 8	0.01	ND	0.033	ND		0.12	0.0012	ND	ND	
R5. 8. 25					ND					
R5. 8. 30					ND					
R5. 9. 3	0.01	ND	0.005		ND	0.13	ND	ND	ND	
R5. 9. 12					ND					
R5. 9. 19					ND					
R5. 9. 26					ND					
R5. 10. 8					ND					
R5. 10. 12	0.02	ND	0.029		ND	ND	0.0010	ND	ND	
R5. 10. 20					ND					
R5. 10. 24					ND					
R5. 11. 3					ND					
R5. 11. 9	0.02	ND	0.014		ND	0.32	0.0011	ND	0.010	
R5. 11. 14					ND					
R5. 11. 22					ND					
R5. 11. 28					ND					
R5. 12. 5	0.02	ND	0.009		ND	0.06	0.0006	ND	ND	
R5. 12. 15					ND					
R5. 12. 20					ND					
R6. 1. 10					ND					
R6. 1. 18	0.02	ND	0.007		ND	0.08	0.0007	ND	ND	
R6. 1. 31					ND					
R6. 2. 9	0.02	ND	0.012		ND	0.05	ND	ND	ND	
R6. 2. 15					ND					
R6. 3. 15	0.02	ND	0.019		ND	0.58	0.0009	ND	0.007	